



UPSTATE

MEDICAL UNIVERSITY

Materials Management


Clinical Supplies, Materials and Logistics



ETCO2 PHILLIPS REPLACED WITH SALTER LABS LAWSON # 60086

Current
Compound Benzoin Tincture

SUB
Compound Benzoin Tincture



Directions for Use

1. Clean and dry skin for at least 2 ½ inches (6.4 cm) around the wound making sure to remove skin oil or exudate. For best results, use a sponge moistened with alcohol or saline followed by a firm wipe with a dry sponge.
2. Using aseptic technique, peel back lid and deliver vial.
3. Gently squeeze vial to crush inner ampoule, thereby releasing compound benzoin tincture. Do not roll crushed ampoule between fingers as this may cause glass to penetrate outer vial.
4. Carefully apply a thin film of compound benzoin tincture to the skin on each side of the wound with strokes parallel to the wound. Begin at the wound edges and work outward. Avoid contact with the wound.
5. Allow tincture to dry before applying surgical tapes and/or adhesive wound closure bandages (approximately 30-45 seconds).