



UPSTATE

MEDICAL UNIVERSITY

Materials Management

Clinical Supplies, Materials and Logistics



VITAL AF 1.2 Cal on Backorder (Lawson #: 17021)

Please see below for supplement replacements:

Vital AF 1.2 Cal Pour Into Pump Pouch

Vital AF 1.2
On Backorder

Replacement

Vital AF 1.2 Cal

Must use Pouch
LABEL
Lawson # 54550



Please contact the dietitian
via VOCERA or call
4-4987(UH) or 5440 (CH)

Hang Time: Sterile, commercially available, decanted formulas (cans/carton)- 8 hours. Allow the bag to completely empty. Prior to refilling, clean top of can/carton with alcohol before opening, rinse the bag with sterile water and re-prime the bag and tubing with fresh formula before connecting to the feeding tube.

We will continue to source our usual supplies as they are available, and return to normal as the supply chain stabilizes. Our suppliers share this disruption may last for several months. Thank you for your understanding and patience.